



NTC8382 and NTC8385 Expression Vectors Instruction Manual

Catalog Numbers

NTC-DV8382
NTC-DV8385-EGFP
NTC-DV8382VA1
NTC-DV8385VA1-EGFP
NTC- DV8382-LV
NTC- DV8385-LV
NTC- DV8382VA1-LV
NTC- DV8385VA1-LV

Version 3

August 2014

Nature Technology Corporation.

4701 Innovation Drive Suite 103 Lincoln Nebraska, 68521

Telephone: (402) 323-6289

Fax: (402) 323-6292

Email: hodgson@natx.com

Website: www.natx.com

General Information

Contents: 20 µg each of plasmid vector shipped in 1x TE buffer.

Storage: Plasmids should be stored at -20°C

Table of Contents

Contents	ii
Storage	ii
Introduction.....	1
Antibiotic-free Selection.....	1
Chimeric Promoter.....	2
Cloning.....	3
Vector Construction	4
Vector Features	5
Intracellular Targeting	5
Plasmid Expression Vector kit.....	6
Linearized Vector kit	6
Accessory Products	7
Cloning into NTC8382 and NTC8385 vectors	8
Transformation and bacterial propagation overview	10
Transformation Procedure	10
Notes	11
SOC Medium	11
References.....	12
Patent and Licensing Information.....	13
Vector Map, Features, and Sequences	
NTC8382.....	14
NTC8385-EGFP	16
NTC8382-VA1	18
NTC8385-VA1-EGFP	20

NTC8382 and NTC8385 vector family

Introduction

Nature Technology Corporations (NTC's) NTC8382 and NTC8385 plasmids are modified antibiotic-free versions of the corresponding kanamycin resistant NTC7482 and NTC7485 vectors.

These plasmids were specifically designed as safe minimalized antibiotic-free selection vectors for the expression of recombinant proteins in mammalian cells. This may be for protein production, or induction of neutralizing immune responses by genetic immunization. The vectors combine minimal prokaryotic sequences including an antibiotic-free sucrose selectable marker. The vectors also contain a novel chimeric promoter that directs superior mammalian cell expression (Luke *et al.* 2009).

The vectors are available in two transgene-targeting versions. NTC8385 expresses encoded protein without additional sequences. NTC8382 targets encoded protein into the secretory pathway using an optimized tissue plasminogen activator (TPA) signal peptide.

The vectors were designed to be responsive to Food and Drug Administration (FDA) regulatory guidance's regarding DNA Vaccine vector composition (FDA 1996, FDA 2007; reviewed in Williams *et al.* 2009a). All sequences that were not essential for *Escherichia coli* plasmid replication or mammalian cell expression of the target gene were eliminated. Synthetic eukaryotic mRNA leader and polyA sequences were utilized in the vector design to limit DNA sequence homology with the human genome to reduce the possibility of chromosomal integration.


Target gene expression is driven from an optimized chimeric promoter-intron (CMV-HTLV-I R synthetic intron). The chimeric CMV promoter achieves significantly higher expression levels than traditional human cytomegalovirus (CMV) promoter based vectors (Luke *et al.* 2009). The vectors encode a consensus Kozak translation initiation sequence and ATG start codon.

Antibiotic-free Selection

Antibiotic-resistance markers, typically kanamycin resistance (KanR), allow selective retention of plasmid DNA during bacterial fermentation and are the most commonly utilized selectable markers.

To ensure safety, however, regulatory agencies recommend elimination of antibiotic-resistance markers from therapeutic and vaccine plasmid DNA vectors. The presence of an antibiotic resistance gene in the plasmid backbone is considered undesirable by regulatory agencies, due to the potential transfer of antibiotic resistance to endogenous microbial flora and the potential activation and transcription of the genes from mammalian promoters after cellular incorporation into the genome (Reviewed in Williams *et al.* 2009a).

The NTC838 series antibiotic-free DNA vaccines (e.g. NTC8385-EGFP; NTC8385-EGFP-VA1; NTC8382; NTC8382 -VA1) incorporate and express a 150 bp RNA-OUT antisense RNA. RNA-OUT represses expression of a chromosomal counter-selectable marker (*SacB*) (**Fig. 1**; Luke *et al.*, 2009). *SacB* encodes a levansucrase, which is toxic in the presence of sucrose. Plasmid selection is achieved in sucrose-containing media.



Levansucrase (*SacB*) conditional-lethal gene on bacterial chromosome

Figure 1: Antibiotic-free selection

NTC838 series sucrose selectable DNA vaccine vectors combine antibiotic-free selection with highly productive HyperGRO fermentation manufacturing (>0.8 g/L plasmid DNA yields), without compromising *in vivo* expression of encoded proteins (Luke *et al.*, 2009).

NTC8382-VA1 and NTC8385-VA1 are derivatives that incorporate the adenoviral serotype 5 VA RNAI (VA1) transient expression enhancer. This modification further improves eukaryotic expression without affecting *Escherichia coli* production yields.

Chimeric promoter

The chimeric promoter is composed of:

- 1) CMV promoter and start of exon 1;
- 2) A HTLV-I R sequence which contains the 5' splice acceptor site;
- 3) A synthetic 3' acceptor site based on the rabbit β globin intron;
- 4) Exon 2 splicing enhancer comprised of serine-arginine rich (SR) protein binding site (3 copies of GAAGAAGAC) to improve RNA export (Lavigneur *et al.* 1993);
- 5) Exon 2 kozak sequence upstream of the start codon for the gene of interest.

Incorporation of the HTLV-I R region downstream of the CMV promoter has been demonstrated to improve expression and cellular immune responses to HIV DNA vaccines in mice and nonhuman primates compared to CMV promoter based vectors (Barouch *et al.* 2005). Expression levels of EGFP from plasmids incorporating the NTC chimeric promoter are also dramatically increased versus CMV promoter containing plasmids (Luke *et al.* 2009; **Fig. 2**).

The HTLV-I R region in NTC838 series vectors (incorporated as part of Exon 1 and Intron 1) dramatically increases mRNA translation efficiency, but not overall mRNA levels, after transient transfection. A similar mRNA translation efficiency increase also occurs with plasmid vectors incorporating and expressing the protein kinase R (PKR)-inhibiting Adenoviral VA RNAI (VA1; NTC8385-EGFP-VA1; NTC8382 -VA1). Strikingly, the HTLV-I R region and VA1 RNA do not increase transgene expression (or mRNA translation efficiency) from plasmid DNA after forced genomic integration.




Figure 2: NTC8385-VA1 expression compared to gWIZ-EGFP and pVAX1-EGFP in HEK293 cells

Cloning

For precise cloning, genes are copied by PCR amplification from clones or genomic DNA using primers with *SalI* (5' end) and *BglII* (3' end) sites. Alternatively, genes are synthesized chemically to be compatible with the *SalI* / *BglII* cloning sites. Design criteria for gene synthesis have been developed by NTC (Williams *et al.* 2009).

In summary, the NTC8382 and NTC8385 vectors offer the following advantages:

- Highest mammalian cell expression using a chimeric CMV promoter-synthetic intron
- Antibiotic-free selection in *Escherichia coli* host
- Optional N-terminal TPA secretion tag for protein export (NTC8382)
- Optional VA1 RNA-based transient expression enhancer
- Simultaneous cloning into all vectors through use of compatible precision cloning cassettes
- Small vectors for more efficient transfection
- Compliance with regulatory guidance (*i.e.* Reduced size, reduced human genome homology, elimination of antibiotic resistance marker)

NTC8382 and NTC8385 vector construction

NTC used Gene Self-Assembly (GENSA) technology to create a series of validated modular elements for vector assembly. Each modular element was assigned a position in a circular array, providing a promoter, 5-leader/splice-site, target gene or high throughput cloning site, terminators, and prokaryotic origin/selection/terminator sites. Individual modules were assigned a position and directionality by means of 4bp unique, non-palindromic address tags. Modules were constructed using a novel cloning process that generates unique, non-palindromic address tags. This allowed precise minimal vector design to eliminate all extraneous sequences.

The pDNAVACCultra vectors (Williams *et al.* 2006) were assembled from GENSA modules consisting of products representing:

- Optimized inducible high copy number pUC prokaryotic replication origin;
- Prokaryotic selectable marker gene (kanamycin);
- Eukaryotic enhancer-promoter (CMV)
- Optimized synthetic eukaryotic untranslated leader-intron-translational initiation sequence (Kozak sequence) cassette
- Gene leader cassette containing the cloning site;
- Optimized synthetic eukaryotic polyA.

NTC8382 and NTC8385 vector features

Feature	Function
Optimized promoter upstream region (UP)	Sequence optimized to maximize CMV promoter expression
Optimized human cytomegalovirus (CMV) immediate-early enhancer promoter	High-level mammalian cell expression <i>in vitro</i> and <i>in vivo</i>
HTLV-I R	Increase mammalian cell expression
Synthetic Intron	Increase mRNA nuclear export
SR-protein binding sites (3x)	Increase mRNA nuclear export
Precision cloning cassette	Clone target gene into the vector with native or secreted intracellular trafficking
Synthetic eukaryotic polyadenylation signal- polyA	Effective mRNA transcriptional termination and polyadenylation
<i>trpA</i> prokaryotic terminator	Protection of replication origin from insert initiated transcription
Adenoviral Serotype 5 VA RNAI (VA1) (Optional)	Increased mammalian cell expression
pUC replication origin	High copy number plasmid production in <i>Escherichia coli</i> cells
Sucrose selection marker (RNA-OUT)	Plasmid selection in <i>Escherichia coli</i> cells*

* Selection in engineering host strain DH5 α att λ ::P_{5/6} 6/6-RNA-IN- SacB, catR

NTC8382 and NTC8385: Intracellular targeting

Vector	Protein Destination	Targeting Tag
NTC8382	secreted	human tissue plasminogen activator (TPA)
NTC8385	native	ATG*

* Native vectors express the target gene from a vector encoded ATG start codon immediately downstream of an optimized kozak sequence (*SalI* site). Protein targeting (*e.g.* cytoplasmic or nuclear) will be determined by protein-intrinsic factors

The TPA targeting peptide has been demonstrated to efficiently target heterologous proteins for secretion (Zhongming *et al.* 1999).

Plasmid Expression Vector kit‡

Vector	Targeting	Quantity	Catalog Number	Price
NTC8382	secreted	20 µg	NTC- DV8382	\$420.00
NTC8385-EGFP*	cytoplasmic EGFP	20 µg	NTC- DV8385-EGFP	\$420.00
NTC8382-VA1	secreted	20 µg	NTC- DV8382VA1	\$420.00
NTC8385-VA1-EGFP*	cytoplasmic EGFP	20 µg	NTC- DV8385VA1-EGFP	\$420.00

* The NTC8385-EGFP control plasmid is used as a transfection control for expression in a cell line of interest, and as the NTC8385 backbone for cloning new genes

‡ Kit contains plasmid and DH5α att_λ::P_{5/6 6/6}-RNA-IN- SacB, catR host strain glycerol stock

Linearized Vector kit‡

Vector	Targeting	Quantity	Catalog Number	Price
NTC8382	secreted	1 µg	NTC- DV8382-LV	\$444.00
NTC8385	native	1 µg	NTC- DV8385-LV	\$444.00
NTC8382-VA1	secreted	1 µg	NTC- DV8382VA1-LV	\$444.00
NTC8385-VA1	native	1 µg	NTC- DV8385VA1-LV	\$444.00

‡ *SalI* - *BglIII* linearized vector sufficient for 20 cloning reactions and DH5α att_λ::P_{5/6 6/6}-RNA-IN- SacB, catR host strain glycerol stock

RAPID-VACCtm CONTRACT CLONING is an optional service whereby NTC clones investigator-specified genes into the desired expression vectors, providing sequence-validated clones, or purified plasmid DNA ready for animal testing (see www.natx.com).

RAPID-VACCtm

Seamless cloning gene of interest into any vector	\$400.00
1mg endotoxin-free plasmid preparations, each	\$335.00

NTC offers the following products for use the NTC8382 and NTC8385 vectors

Accessory Products

Product	Description	Quantity	Catalog Number	Price
pVAC5'	Forward sequencing primer	500 pmol (10 pmol/μL)	NTC-DVU-SP1	\$30.00
pVAC3'	Reverse sequencing primer	500 pmol (10 pmol/μL)	NTC-DVU-SP2	\$30.00
NTC4862 DH5α att _λ ::P _{5/6 6/6} -RNA-IN- SacB, catR electrocompetent cells	Host strain for NTC8382 and NTC8385 selection	0.5 mL	NTC-DVU-CC1	\$199.00
NTC48165 DH5α dcm att _λ ::P _{5/6 6/6} -RNA-IN- SacB, catR electrocompetent cells	dcm- host strain for NTC8382 and NTC8385 selection	0.5 mL	NTC-DVU-CC2	\$199.00
Sucrose selection agar	Plates for NTC8382 and NTC8385 selection	Pouch to make 1 L agar*	NTC-DVU-MD1	\$20
Sucrose selection media	Media for NTC8382 and NTC8385 culture	Pouch to make 1 L media*	NTC-DVU-MD2	\$20

* Reconstitute with distilled water to 1L and autoclave prior to use

Cloning into NTC8382 and NTC8385 vectors

Overview: An example strategy for cloning into the NTC8382 and NTC8385 vectors is outlined below.

GTCGACATG----- Gene of interest----Stop codon -----AGATCT
SalI *BglIII*

For both vectors, the ATG start codon (double underlined) is immediately preceded by a *SalI* site. In NTC8385, the *SalI* site has been demonstrated to be an effective consensus Kozak sequence for translational initiation. In NTC8382, the *SalI* site is downstream in frame with the optimized TPA secretion sequence.

For precise cloning, genes are copied by PCR amplification from clones or genomic DNA using primers with *SalI* (5' end) and *BglIII* (3' end) sites. Alternatively, genes are synthesized chemically to be compatible with the *SalI* / *BglIII* cloning sites. Design criteria for gene synthesis are reviewed in Williams *et al.* 2009a.

For NTC8385, the start codon ATG must immediately follow the *SalI* site (GTCGACATG) or, optionally, a kozak site (*e.g.* GTCGACGCCACCATG). For NTC8382, the ATG is optional but the same reading frame must be retained. For both vectors one or two stop codon (preferably TAA or TGA) must be included prior to the *BglIII* site. A PCR product designed for NTC8385 is compatible with, and can also be cloned into, the NTC8382 vector.

Protocol: Genes are PCR amplified with primers incorporating a *SalI* site into the 5' termini and a *BglIII* site into the 3' termini. Example primers, for cloning a target gene of interest into both vectors are diagramed below.

SalI
GCGCATGGTCGACATG-GENE OF INTEREST 5' END
BglIII
GCAGAAAGATCTTTA-GENE OF INTEREST (REVERSE COMPLEMENT) 3' END

Cleavage of the vectors with *SalI* / *BglIII* generates sticky ends compatible with the cleaved PCR product. The insert is thus directionally and precisely cloned into the vector. The majority of recovered colonies are recombinant, since the generated sticky ends in the parental vector are not compatible.

If the gene of interest encodes internal *Bgl*III sites, *Bam*HI or *Bcl*II restriction enzymes can be utilized in the PCR primers since they create *Bgl*III compatible sticky ends. Note that both enzymes sites will be lost in the resulting construct.

*Xho*I is not recommended as an alternative to *Sal*I for cloning into NTC8385, since the ligated hybrid (GTCGAGATG) has not been established as an effective kozak sequence. This is not an issue for cloning downstream of TPA in NTC8382.

Alternatively, when these restriction sites are present in the gene, the PCR primers may be modified to incorporate *Aar*I (Fermentas, Vilnius Lithuania) type IIS enzyme sites. This enzyme cleaves after +4 and +8 bp downstream of the *Aar*I recognition sequence to create any designed 4 bp 3'-recessed sticky end. *Aar*I primer design to amplify a gene for cloning into NTC8382 and NTC8385 vectors is shown below.

*Aar*I

5' CTCCAGCACCTGCCTATTCGACATG-GENE OF INTEREST 5' END

5' CGTGAGCACCTGCAACGGATCTTA-GENE OF INTEREST (REVERSE COMPLEMENT)
3' END

Recombinant clones can be identified by restriction digestion. *Sal*I and *Bgl*III release the gene insert when using *Aar*I or *Sal*I / *Bgl*III containing primers. *Bgl*III / *Bam*HI / *Bcl*II ligated cohesive termini in the resultant clone will not cleave with either parent restriction enzyme.

The sequence of the gene insert can be determined using the following cloning cassette flanking sequencing primers.

pVAC5': GCTTTTCTGCCAGGTGCTGA
(hybridizes to intron and sequences from 5' end of gene)

pVAC3': GCCAGAAGTCAGATGCTCAA
(hybridizes to polyA region and sequences from 3' end of gene)

Transformation and bacterial propagation overview

RNA-OUT represses expression of a sucrose counter-selectable marker (*SacB*) from the host chromosome. Ligation reactions must therefore be transformed into *SacB* expression strain DH5 α att λ ::P_{5/6 6/6}-RNA-IN- *SacB*, *catR* competent cells.

Cells are plated and propagated on 6% sucrose media (Note: **Plates must not contain NaCl**)

A 50% Sucrose solution should be prepared (wt./vol.) prior to making liquid or solid media. (This should not be autoclaved but rather filter-sterilized with a 0.2 micron filter prior to use)

Solid Media Preparation	Liquid Media Preparation
10g tryptone 5 g yeast extract 15 g agar QS to 880 mL with H ₂ O Autoclave 20 min 121°C, when cooled to ~50°C, add 120 mL 50% Sucrose (wt./vol.), mix, pour plates.	10g tryptone 5 g yeast extract QS to 880 mL with H ₂ O Autoclave 20 min 121°C When cooled to ~50°C, add 120 mL 50% Sucrose (wt./vol.), allow to cool before use.

Transformation Procedure

1. Ligation reactions must be diluted (See note 1) to avoid arcing in the electroporator. Dilute reactions 5-fold with sterile water and add 1ul to a chilled microcentrifuge tube. Prechill on ice one 0.1 cm electroporation cuvette for each sample to be electroporated.
2. Thaw NTC electrocompetent cells on wet ice.
3. When cells are thawed, mix cells by tapping gently. Add 20 μ l of cells to each chilled microcentrifuge tube containing your DNA to be transfected.
4. Refreeze any unused cells. Although the cells are refreezable, subsequent freeze-thaw cycles will decrease transformation efficiency.
5. Pipette the cell/DNA mixture into a prechilled 0.1 cm cuvette and electroporate. If you are using the BTX®ECM®630 or BioRad GenePulser®II electroporator, we recommend using the following electroporator conditions: 2.0 kV, 200 Ω , 25 μ F. (See note 2)
6. To the cells in the cuvette, quickly add 0.5-1 ml of room temperature S.O.C. medium and transfer the solution to a 15 ml snap-cap tube (e.g. Falcon™ tube).
7. Shake at 225 rpm (30°C) for 1 hour. (This step may be omitted if transforming supercoiled plasmid DNA)
8. Spread 50-100 μ l on Sucrose selection agar (or dilute experimental reactions as necessary and spread 50-100 μ l on selective plates)
9. Incubate plates overnight at 30°C.

Notes

1. Transformation efficiencies will be 10- to 100-fold lower for ligation mixtures than for an intact control plasmid. **Salts and buffers severely inhibit electroporation.** Ligation reactions can be diluted 5-fold with sterile water, and 1 μ l added to 20 μ l of cells. Adding undiluted ligation mixtures or too high a volume of DNA decreases transformation efficiency and increases the risk of arcing.
2. If you are using an electroporator other than a BTX®ECM®630 or BioRad GenePulser®II electroporator, you may need to vary the settings to achieve optimal transformation efficiency.

SOC Medium

1. Add the following to 900ml of distilled H₂O
 - 20g Bacto Tryptone
 - 5g Bacto Yeast Extract
 - 2ml of 5M NaCl.
 - 2.5ml of 1M KCl.
 - 10ml of 1M MgCl₂
 - 10ml of 1M MgSO₄
 - 20ml of 1M glucose
2. Adjust to 1L with distilled H₂O
3. Sterilize by autoclaving

References

- Barouch DH, Yang ZY, Kong WP, Koriath-Schmitz B, Sumida SM, Truitt DM, Kishko MG, Arthur JC, Miura A, Mascola JR, Letvin NL, and Nabel GJ. (2005) A human T-cell leukemia virus type 1 regulatory element enhances the immunogenicity of human immunodeficiency virus type 1 DNA vaccines in mice and nonhuman primates. *J. Virol.* 79: 8828-8834
- FDA. (1996) Points to consider on plasmid DNA vaccines for preventive infectious disease indications. US Food and Drug Administration
- FDA. (2007) Guidance for industry: considerations for plasmid DNA vaccines for infectious disease indications. US Food and Drug Administration
- Lavigueur A, La Branche H, Kornblihtt AR, and Chabot B. (1993) A splicing enhancer in the human fibronectin alternate ED1 exon interacts with SR proteins and stimulates U2 snRNP binding. *Genes Develop.* 7: 2405-2417
- Luke J, Carnes AE, Hodgson CP, and Williams JA. (2009) Improved antibiotic-free DNA vaccine vectors utilizing a novel RNA based plasmid selection system. *Vaccine* 27: 6454-6559
- Williams JA, Carnes AE, and Hodgson CP. (2009a) Plasmid DNA vector design; impact on efficacy, safety and upstream production. *Biotechnology Advances* 27:353-370
- Williams JA, Luke J, Langtry S, Anderson S, Hodgson CP, and Carnes AE. (2009b) Generic plasmid DNA production platform incorporating low metabolic burden seed-stock and fed-batch fermentation processes. *Biotechnol Bioeng* 103:1129-1143
- Williams JA, Luke J, Johnson L, and Hodgson CP. (2006) pDNAVACCultra vector family: high throughput intracellular targeting DNA vaccine plasmids. *Vaccine* 24:4671-4676
- Williams JA. (2014) Improving DNA vaccine performance through vector Design. *Curr Gene Ther* 14:170-189
- Zhongming L, Howard A, Kelley C, Delogu G, Collins F, and Morris S. (1999) Immunogenicity of DNA vaccines expressing tuberculosis proteins fused to tissue plasminogen activator signal sequences. *Infect Immun.* 67: 4780-4786

Patent and Licensing information

Limited License

Nature Technology Corporation (NTC) grants the end user (purchaser) of the NTC8382 and NTC8385 expression vectors and NTC4862, NTC48165, NTC5402, and NTC54208 sucrose-selection host strains a nontransferable, non-exclusive license to use the plasmids and host strains for non-commercial research purposes only. These vectors and host strains are intended for research use only by the purchaser.

The purchaser cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for commercial purposes.

Separate licenses are available from NTC for the express purpose of non-research use or applications of the NTC8382 and NTC8385 vectors and NTC4862, NTC48165, NTC5402, and NTC54208 host strains.

Product Use Limitations

The NTC8382 and NTC8385 vectors and NTC4862, NTC48165, NTC5402, and NTC54208 host strains are sold for research use only. They are not to be used for human diagnostic or included/used in any drug intended for human use.

Patent Information


The CMV promoter is covered under U.S. Patents 5,168,062 and 5,385,839 and its use is permitted for research purposes only. Any other use of the CMV promoter requires a license from the University of Iowa Research Foundation, 214 Technology Innovation Center, Iowa City, IA 52242. The polymerase chain reaction (PCR) process is covered by patents owned by Roche and requires a license for use.

The vectors and cell lines herein are covered by World Patent Application WO2008153733, JA Williams, Vectors and methods for genetic immunization (2008), and by US, European, and Australian patents: US 2010/0303859; EP2333091; and AU 2008262478, respectively.

NTC makes no representations that the use of the NTC8382 and NTC8385 vectors and/or NTC4862, NTC48165, NTC5402, and NTC54208 host strains will not infringe any patent, copyright, trademark, or other proprietary rights.

For more information, please contact: Clague Hodgson
Nature Technology Corporation, 4701 Innovation Drive, Suite 103, Lincoln Nebraska, 68521
Telephone: (402) 323-6289, **Fax:** (402) 323-6292, **Email:** natx@natx.com.

NTC8382




Vector Features

trpA prokaryotic terminator: 2897-28
pUC replication origin: 35-1049
Sucrose selection marker (RNA-OUT): 1066-1210
Promoter Upstream Region (UP): 1225-1518
CMV enhancer: 1519-2073
CMV promoter: 2074-2193
Untranslated leader (exon 1): 2194-2380
HTLV-I R: 2265-2490
Synthetic Rabbit β -globin-based 3' intron: 2499-2605
Intron: 2381-2605
Exon 2 (SR-protein binding sites-Kozak):2606-2652
TPA N-terminal targeting tag: 2653-2721
Sall-BglIII cloning cassette: 2734-2751
PolyA: 2746-2886

ccgcctaagagcgggcttttttcttaggccttctccgcttccctgctcactgactcgtcgcgctcggctgctcggcga
gcggtatcagctcactcaaaggcggtaatacggttatccacagaatcaggggataacgcaggaaagaacatgtgagcaaaag
gccagcaaaaggccaggaaccgtaaaaaaggccgcgttgctggcgttttccataggtccgccccctgacgagcatcaca
aatcgacgctcaagtcagaggtggcgaacccgacaggactataaaagataaccaggcggttccccctggaagctccctcgtgcg
ctctctgttccgacctgccgcttaccggatacctgtccgctttctccctcgggaagcgtggcgctttctatagctcacgctgt
aggtatctcagttcgggtgtaggtcgtcgtccaagctgggctgtgtgcacgaacccccgttcagcccagccgctgcgccttat
ccggaactatcgtctgagccaacccgtaagacacgacttatgccactggcagcagccactggtaacaggattagcagag
cgaggtatgtaggcgggtgtacagagtcttgaagtgggtggcctaactacggctacactagaagaacagtatttggtatctgcgct
ctgctgaagccagttacctcggaaaaagattgtagctcttgatccggcaacaaccaccgctgtagcgggtggtttttgtt
tgcaagcagcagattacgcgcagaaaaaaaggatctcaagaagatcctttgatctttctacggggtctgacgctcagtggaacg
aaaactcacgtaagggtatttgatgagattatcaaaaaggatctcactagatcctttaaataaaatgaagtttaaatcaa
tctaaagtatatatgagtaactgtgctgacagttaccaatgcttaacagtgaggcacctatctcagcgtatctgtctattcgttcat
ccatagttgctgactcctgcaaacacggttggtgtagaattggtaaagagagtcgtgtaaaatcagagttcgacatctgttgtc
tgattattgattttggcgaaaccattgatcatatgacaagatgtgatctaccttaactaatgatttgataaaaatcattaggtacc
cggctctagatggcatgacattaacctataaaaataggcgtatcacgaggcccttcgtctcgcgcgtttcgggtgatgacgggtgaa
aacctctgacacatgacgctcccggagacggtcacagctgtctgtaagcggatgccgggagcagacaagcccgtcagggcg
cgtcagcgggtgttggcgggtgtcggggctggcttaactatcggcgcacagagcagattgtactgagagtgccacatgctgggt
gtgaaataccgcacagatgcgtaaggagaaaataccgcatcagattggctattggccattgcatacgttgatccatataat
gtacatttatattgctcatgtccaacattaccgccatgtgacattgattattgactagtattaatagtaataacacggggcatta
gttcatagccatataatggagttccgcttacataactacggtaaatggcccgcctggctgaccgccaacgacccccgccat
tgacgtcaataatgacgtatgtccatagtaacgcaatagggactttccattgacgtcaatgggtggagtatttacggtaaacg
cccacttggcagfacatcaagtgtatcatatgcaagtagccccctattgacgtcaatgacggtaaatggcccgcctggcattat
gcccagfacatgaccttatgggactttcctacttggcagfacatctacgtattagtcacgctattaccatgggtgatgcggtttggca
gtacatcaatgggcgtggatagcgggttgactcacggggattccaagtctccaccctgacgtcaatgggagtttgtttggca
ccaaaatcaacgggactttccaaaatgtcgtacaactccgcccattgacgcaaatgggcggtaggcgtgtacgggtgggaggt
ctatataagcagagctcgttttagtaaccgtcagatcgctggagacgccatccacgctgtttgacctcatagaagacaccgg
gaccgatccagcctccgcggctcgcactctctctcaccgcgccgcccctacctgaggccgcatccacgccggtgagtc
gcgttctgccgctcccgcctgtggtgcctctgaactcgcgtccgccgtctaggttaagtttaagctcaggtcagaccggcct
ttgtccggcgtcccttggagcctacctagactcagccggctctccacgctttgctgacctgcttgcactctagttctctcgt
taactfaatgagacagatagaaactggtctttagaacaacagagtagtcgctgcttttctgccaggtgctgacttctctcccctggg
ctttttcttttctcaggtgaaaagaagaagacgaagaagacgaagaagacaaagccgccaccatggatgcaatgaagagag
ggctctgctgtgtgctgctgtgtggagcagctcttctgttccgcccagcggtagccgatccgtcaggggggagatcttttcc
ctctgcaaaaattatggggacatcatgaagccccttgagcatctgacttctggctaataaaggaaatttatttcattgcaatagtg
gttgaattttgtgtctctcactcggaggacataagggcggccgctagc

NTC8385- EGFP




Vector Features

trpA prokaryotic terminator: 3530-28
 pUC replication origin: 35-1049
 Sucrose selection marker (RNA-OUT): 1066-1210
 Promoter Upstream Region (UP): 1225-1518
 CMV enhancer: 1519-2073
 CMV promoter: 2074-2193
 Untranslated leader (exon 1): 2194-2380
 HTLV-I R: 2265-2490
 Synthetic Rabbit β -globin-based 3' intron: 2499-2605
 Intron: 2381-2605
 Exon 2 (SR-protein binding sites-Kozak): 2606-2654
 EGFP: 2655-3378
 PolyA: 3379-3519

ccgcctaagagcgggcttttttcttaggccttctccgcttccctgctcactgactcgtcgcgctcggctgctcggcga
gcggfatcagctcactcaaaggcggtaatacggttatccacagaatcaggggataacgcaggaaagaacatgtgagcaaaag
gccagcaaaaggccaggaaccgtaaaaaaggccgcgttgctggcgttttccataggtccgccccctgacgagcatcaca
aatcgacgctcaagtcagaggtggcgaaccgacaggactataaaagataaccaggcggttccccctggaagctccctcgtc
ctctctgttccgacctgccgcttaccggatacctgtccgcttctccctcgggaagcgtggcgcttctcagctcacgctgt
aggtatctcagttcgggtgtaggtcgtcgtccaagctgggctgtgtgcagcaacccccgttcagcccagccgctgcgctt
ccggaactatcgtctgagccaaccggtaagacacgacttatgccactggcagcagccactggaacaggattagcagag
cgaggtatgtaggcgggtgtacagagtcttgaagtgggtggcctaactacggctacactagaagaacagtatttggtatctgcgct
ctgctgaagccagttacctcggaaaaagattgtagctctgatccggcaacaaccaccgctgtagcgggtggtttttgtt
tgcaagcagcagattacgcgcagaaaaaaggatctcaagaagatcctttgatctttctacggggtctgacgctcagtggaac
aaaactcacgtaagggtatttgctcatgagattataaaaaggatctcactagatcctttaaataaaatgaagtttaaatcaa
tctaaagtatatatgagtaactgtgctgacagttaccaatgctaatcagtgaggcacctatctcagcgtatctgtctatttcgct
ccatagttgctgactcctgcaaacacggtgtgtagaattgtaaaagagctgtgtaaaatcaggttcgacatctgttgtc
tgattattgattttggcgaaaccattgatcatatgacaagatgtgtatctaccttaactaatgatttgataaaaatcattaggtacc
cggctctagatggcatgacattaacctataaaaataggcgtatcacgaggcccttcgtctcgcgctttcgggtgatgacgggtgaa
aacctctgacacatgacgctcccggagacggtcacagctgtctgtaagcggatgccgggagcagacaagcccgtcagggcg
cgtcagcgggtgttggcgggtgtcggggctggcttaactatcggcgcacagagcagattgtactgagagtgccacatgctggt
gtgaaataccgcacagatgcgtaaggagaaaataccgcatcagattggctattggccattgcatacgtgtatccatataat
gtacatttatattgctcatgtccaacattaccgccatgtgacattgattattgactagtattaatagtaataacacggggtcatta
gttcatagccatataatggagtccgcgttacataactacggtaaatggcccgcctggctgaccgccaacgacccccgccat
tgacgtcaataatgacgtatgtccatagtaacgcaatagggacttccattgacgtcaatgggtggagtatttacggtaaacg
cccacttggcagfacatcaagtgtatcatatgcaagtagccccctattgacgtcaatgacggtaaatggcccgcctggcattat
gcccagfacatgaccttatgggacttccacttggcagfacatctacgtattagctcgtcattaccatggtgatgcggttttggca
gtacatcaatggcggtgatagcgggttactcaggggattccaagtctccaccctgacgtcaatgggagtttgttttggca
ccaaaatcaacgggacttccaaaatgtcgtacaactccgcccattgacgcaaatggcggtgtaggcgtgtacggtgggaggt
ctatataagcagagctcgttagtgaaccgtcagatcgctggagacgccatccacgctgtttgacctcatagaagacaccgg
gaccgatccagcctccgcggctcgcactctctctcaccgcccgcctaccctgagggcccatccacgcccgttgatc
gcgttctgccgctcccgcctgtggtgctcctgaactcgcgtccgcttaggtaagttaaagctcaggtcagaccgggct
ttgtccggcgtcccttggagcctacactagactcagccgctcctccagctttgctgacctgcttcaactctagtctctcgt
taactaatgagacagatgaaaactggtctgtagaacagagtagtcgctgcttttctgccaggtgctgacttctccccctggg
ctttttcttttctcaggttgaagaagaagacgaagaagacgaagaagacaaccgctcgtcagatggtgagcaaggcga
ggagctgtcaccggggtggtgccatcctggtcagctggacggcgacgtaaacggccacaagttcagcgtgtccggcgag
ggcgagggcgatgccactacggcaagctgacctgaagtctctgcaccaccggcaagctgcccgtgcccggcccacc
tcgtgaccacctgacctacggcgtgagcttccagccgctaccccgaccacatgaagcagcagacttcttaagtccgcca
tgcccgaaggctacgtccaggagcgcacctctctcaaggacgacggcaactacaagaccgcccggaggtgaagttcga
ggcgacaccctggtgaaccgcatcagctgaagggcctcactcaaggaggacggcaacatcctggggcacaagctgga
gtacaactacaacagccacaacgtctatatcatggccgacaagcagaagaacggcatcaaggtgaactcaagatccgccaca
acatcgaggacggcagcgtgcagctcggcaccactaccagcagaacccccatcggcgacggccccgtgctgctgccc
acaaccactacctgacaccagtcgcccagcaagacccccagagaagcgcgacacatggtcctgctgaggtcgtg
accgcccccggatcactctcggcatggcagagctgtacaagtaaccgagatcttttccctctgcaaaaattatggggacat
catgaagccccctgagcatctgactctggctaataaaggaaatttttcttcaatgcaatagtggttgaatttttgtgctctcactc
ggaaggacataagggcgccgctagc

NTC8382-VA1




Vector Features

trpA prokaryotic terminator: 3091-28
 VA1 RNA transient expression enhancer: 32-225
 pUC replication origin: 229-1243
 Sucrose selection marker (RNA-OUT): 1260-1404
 Promoter Upstream Region (UP): 1419-1712
 CMV enhancer: 1713-2267
 CMV promoter: 2268-2387
 Untranslated leader (exon 1): 2388-2574
 HTLV-I R: 2459-2684
 Synthetic Rabbit β -globin-based 3' intron: 2693-2799
 Intron: 2575-2799
 Exon 2 (SR-protein binding sites-Kozak): 2800-2846
 TPA N-terminal targeting tag: 2847-2915
Sall-BglIII cloning cassette: 2928-2945
 PolyA: 2940-3080

ccgcctaagagcgggcttttttcttagggcgcaaaaggagagcctgtaagcgggactctccgtggtctggtgataaattc
caagggatcatggcggacgaccggggtcagccccgatccggccgtccgccgtgatccatgcggtaccgcccgcgtgc
gaaccaggtgtgcacgtcagacaacgggggagtgctcctttggcttctcccttctccgcttctcgtcactgactcgt
gcgctcgggtcgtcgggtcggcgagcgggtatcagctactcaagggcggtaatacgggtatccacagaatcaggggataacg
caggaaaacatgtgagcaaaaggccagcaaaaggccaggaaccgtaaaaaggccgcgttgctggcgttttccataggct
ccgccccctgacgagcatcaaaaaatcgacgctcaagtcagaggtggcgaaccggacagactataaagataaccaggc
gtttcccccctggaagctccctcgtgcgtctcctgtccgacctgccgttaccggatacctgtccgccttctccctcgggaag
cgtggcgtttctcatagctcacgctgtaggtatctcagttcgggtgtaggtcgttcgctcaagctgggctgtgtgcacgaacccc
ccgttcagcccagcgtgcgcttatccggttaactatcgtcttgagtcacaccggtaagacacgacttatcgcactggcagc
agccactggtaacaggattagcagagcgggtatgtagggcgtgtacagagttctgaagtgggtggcctaactacggctacac
tagaagaacagtatgtgtagctcgtcgtcgaagccagttacctcggaaaaagagttgtagctcttgatccggcaaaaa
ccaccgctgtagcgggtgtttttgttgaagcagcagattacgcgcagaaaaaggatctcaagaagatcctttgatctttc
tacggggtctgacgctcagtggaacgaaaactcacgttaagggattttggtcatgagattatcaaaaaggatcttcactagatcct
tttaataaaaaatgaagtttaaatcaatctaaagtatatatgagtaaacctggctctgacagttaccaatgcttaacagtgaggcacc
tatctcagcgtatctgtctatttcgttcatccatagttgcctgactcctgcaaacacggttggtgtagaattggtaaagagagtcgtgta
aaatcgcagttcgacatcttgtgtctgattattgattttggcgaaaccattgatcatatgacaagatgtgatctaccttaacttaa
tgattttgataaaaaatcattaggtacccggctctagatggcatgacattaacctataaaaaataggcgtatcacagggcccttctcgtc
tcgctcgtttcgggtgatgacgggtgaaaacctctgacacatgcagctccggagacggtcacagctgtctgtaagcggatgccg
ggagcagacaagcccgtcagggcgcgtcagcgggtgtggcgggtgtcggggctggcttaactatcgggcatcagagcagat
tgactgagagtgaccatagcgggtgaaataccgcagatgcgtaaggagaaaaataccgcatcagattggctattggccat
tgcatacgttgatccatataatgtacattatattggctcatgccaacattaccgcatgttgacattgattattgactagtatt
aatagtaatacattacggggcattagttcatagccatataaggagttccgcgttacataacttacggtaaattggccgctggct
gaccgccaacgacccccgccattgacgtcaataatgacgtatgtccatagtaacgccaatagggactttccattgacgtca
atgggtggagatttacgggtaaacctgccactggcagtagatcaagtgatcatatgccaagtacccccctattgacgtcaatga
cggtaaatggccgctggcattatgccagtagacattatgggacttctacttggcagtagatctacgtattagtcacgct
attaccatggtgatcgggtttggcagtagatcaatgggctgtgtagcgggttgactcacggggatttccaagtctccacccatt
gacgtcaatgggagttgtttggcaccaaaatcaacgggactttccaaaatgtcgtacaactccgccccattgacgcaaatggg
cggtaggcgtgtacgggtgggaggtctatataagcagagctcgttagtgaaccgtcagatgcctggagacccatccacgctg
ttttgacctccatagaaacaccgggaccgatccagcctccgcggctcgcactctccttcacgcgccgcccctacctgag
gcccacatccacgcccgtgagtcgcttctgcccctcccgcctgtggtgcctcctgaactgcgtccgccgtctaggttaagtt
aaagctcaggtcagaccgggctttgtccggcgtcccttggagcctacctagactagccggctcaccgctttgcctgacc
ctgcttctcaactctagttctctgftaacttaatgagacagatagaaactggctttagaaacagagtagtcgctgctttctgc
caggtgctgacttctctcccctgggcttttttcttctcaggtgaaaaagaagacgaagaagacgaagaagacaaaagccgc
caccatggatgcaatgaagagagggctctgctgtgtgctgctgtgtgtagcagctctcgttccagcgggtaccggatcc
gtcagggggggagatcttttccctctgcaaaaaatfatggggacatcatgaagccccttgagcatctgacttctggctaataaa
ggaaattttttcattgcaatagtggttgaatttttgtgtctcactcgggaaggacataagggcgccgctagc

NTC8385-VA1- EGFP



Vector Features

trpA prokaryotic terminator: 3724-28
 VA1 RNA transient expression enhancer: 32-225
 pUC replication origin: 229-1243
 Sucrose selection marker (RNA-OUT): 1260-1404
 Promoter Upstream Region (UP): 1419-1712
 CMV enhancer: 1713-2267
 CMV promoter: 2268-2387
 Untranslated leader (exon 1): 2388-2574
 HTLV-I R: 2459-2684
 Synthetic Rabbit β -globin-based 3' intron: 2693-2799
 Intron: 2575-2799
 Exon 2 (SR-protein binding sites-Kozak): 2800-2848
 EGFP: 2849-3572
 PolyA: 3573-3713

ccgcctaataagcgggcttttttcttagggcgcaaaaggagagcctgtaagcgggactctccgtggtctggtgataaattc
caagggatcatggcggacgaccggggtcagccccgatccggccgtccgccgtgatccatgcggtaccgcccgcgtgc
gaaccaggtgtgcacgtcagacaacgggggagtgctcctttggttcttcccccttctccgcttctcgtcactgactcgt
gcgctcgggtcgtcgggtcggcgagcgggtatcagctactcaaggcgtaatacgggtatccacagaatcaggggataacg
caggaaaacatgtgagcaaaaggccagcaaaaggccaggaaccgtaaaaaggccgcttctggcggttttccataggct
ccgccccctgacgagcatcaaaaaatcagcgtcaagtcagaggtggcgaaccgacaggactataaagataaccaggc
gttccccctggaagctccctcgtgcgtctcctgtccgacctgccgttaccggatacctgtccgccttctccctcgggaag
cgtggcgcttctcatagctcacgctgtaggtatctcagttcgggtgtaggtcgttcgctcaagctgggctgtgtgcacgaaccc
ccgttcagcccagcgtgcgcttatccgtaactatcgtcttgagtcacaaccggtaagacacgacttatcgcactggcagc
agccactggtaacaggattagcagagcgggtatgtaggggtgtacagagttctgaagtgggtggcctaactacggctacac
tagaagaacagtatgtgtagctcgtcgtcgaagccagttacctcggaaaaagagttgtagctcttgatccggcaaaaa
ccaccgctggtagcgggtggttttttgaagcagcagattacgcgcagaaaaaggatcgaagaagatccttggatctttc
tacggggtctgacgctcagtggaacgaaaactcacgttaagggattttggtcatgagattatcaaaaaggatcttcactgatcct
tttaataaaaaatgaagtttaaatcaatctaaagtatatagtaaaactggctctgacagttaccaatgcttaacagtgaggcacc
tatctcagcgtatctgtctatttcgttcatccatagttgcctgactcctgcaaacacggttggtgtagaattggtaaaagagagtcgtgta
aaatcgcagttcgacatcttgtgtctgattattgattttggcgaaaccattgatcatatgacaagatgtgtatctaccttaactta
tgatttgataaaaaatcattaggtacccggctctagatggcatgacattaacctataaaaaataggcgatcacagggcccttctcgtc
tcgcgcttccggtgatgacgggtgaaaactctgacacatgcagctccggagacggtcacagcttctgtaagcggatgccg
ggagcagacaagcccgtcagggcgctcagcgggtgtggcggggtcggggctggcttaactatcgggcatcagagcagat
tgactgagagtgaccatagcgggtgaaataccgcagatgcgtaaggagaaaaataccgcatcagattggctattggccat
tgcatacgtttagatccatataatgtacattatattggctcagtcacaacattaccgcatgttgacattgattattgactagtatt
aatagtaataacacggggcattagttcatagccatataaggagttccgcgttacataacttacggtaaatggcccgctggct
gaccgccaacgacccccgccattgacgtcaataatgacgtatgtccatagtaacgcaatagggactttccattgacgtca
atgggtggagatttacggtaaaactccccactggcagatcacatcaagtgtatcatatgcaagtagccccctattgacgtcaatga
cggtaaatggcccgctggcattatgccagtagacattatgggacttctacttggcagtagatctacgtattagtcacgtc
attaccatggtgatcgggtttggcagtagatcaatgggctggtatagcgggttgactcacggggatttcaagctccacccatt
gacgtcaatgggagttgtttggcaccaaaatcaacgggactttcaaaaatgtcgtacaactccgccccattgacgcaaatggg
cggtaggcgtgtacgggtgggaggtctatataagcagagctcgttagtgaaccgtcagatcgctggagacgccatccacgctg
ttttgacctccatagaaacaccgggaccgatccagcctccgcggctcgcactctccttccacgcgcccgccctacctgag
gcccacatccacgccggtgagtcgcttctcgcctcccgcctggtgctcctctgaactgcgtccgccgtctaggttaagtt
aaagctcaggtcagaccgggcttctcggcgctcccttggagcctacctagactagccggctcaccgcttgcctgacc
ctgcttctcaactctagttctctgtaacttaatgagacagatagaaactggctttagaaacagagtagtcgctgctttctgc
caggtgctgacttctcctccctgggcttttttctcaggtgaaaaagaagagcgaagaagacgaagaagacaaaaccgac
gtcagatggtgagcaaggcgagggagctgtcaccggggtggtgccatcctggtcagctggacggcgacgtaaacggc
cacaagttcagcgtgtccggcgagggcgagggcgatgccacctacggcaagctgacctgaagttcatctgaccaccggca
agctgcccgtgccctggcccacctcgtgaccacctgacctacggcgtgcagtgcttcagccgctaccccgaccacatgaag
cagcacgacttctcaagtccgcatgcccgaaggctacgtccaggagcgcaccatcttctcaaggacgacggcaactacaa
gaccgcccggaggtgaagttcagggcgacacctggtgaaccgcatcagctgaaggcagcacttcaaggaggacgg
caacatcctggggcacaagctggagtacaactacaacagccacaacgtctatatcatggccgacaagcagaagaacggcagc
aaggtgaactcaagatccgccacaacatcaggacggcagcgtgcagctcggcaccactaccagcagaacacccccatcg
gcgacggccccgtgctgctgccgacaaccactacctgagcaccagtcggcctgagcaaaagacccaacgagaagcgcg
atcacatggtcctgctggagttcgtgaccgccggggtatcctcggcatggacgagctgtacaagtaaccgcagatcttttc
cctctgcaaaaaatattggggacatcatgaagcccttgacatctgacttctggctaataaaggaaattttttcattgcaatagtg
tgttgaatttttgtctctcactcgggaaggacataaggcgccgctagc